



Deployable SQUID-based magnetic resonance imaging systems

**Per Magnelind, Andrei Matlashov, Shaun Newman,
Henrik Sandin, Robert Sedillo, Algis Urbaitis,
Petr Volegov, and Michelle Espy**

EUCAS, Sept. 9, 2015

Outline

- Motivation
- Methods
- Modeling
- Results
- Outlook
- Summary



*Single-average slice acquired
inside a shielded room*

Motivation

Magnetic Resonance Imaging (MRI)

- best method for non-invasive imaging of soft tissue anatomy
- saves countless lives each year

Conventional (high-field) MRI

- only in large well-funded medical centers
- is not available in rural settings
- is not deployable to emergency situations or battlefield hospitals


Ultra-low field (ULF) MRI

- pulsed pre-polarization at < 0.3 T
- sensitive Superconducting Quantum Interference Device (SQUID) detection
- greatly relaxed homogeneity
- presence of non-magnetic metal is not an issue
- can be light and made portable




Methods—NMR basics

Conventional NMR:




ULF-NMR: non-adiabatic switching of B_p




Methods—NMR basics

ULF-NMR: adiabatic switching of B_p




ULF-NMR: non-adiabatic switching of B_p



Methods—Gradients

Frequency encoding



Phase encoding





Image (spin density) obtained from 3-D inverse Fourier Transform

Methods—Pulse-sequence



3D Fourier spin-echo imaging sequence.

B_p : pre-polarization field
 ≤ 100 mT, linearly ramped over 70 ms, then
 ramped to zero adiabatically during 20 ms.

B_m : measurement field
 $190\text{--}200$ μT (8.2–8.6 kHz).



G_x (readout): ≤ 248 $\mu\text{T/m}$ ≤ 11 Hz/mm

G_z (phase, in-plane): ≤ 222 $\mu\text{T/m}$ ≤ 9 Hz/mm

G_y (phase, depth): ≤ 47.6 $\mu\text{T/m}$ ≤ 2 Hz/mm


Implemented in a custom LabView VI with a
 time-slot matrix interface.

Methods—Unshielded System Hardware




- 3 pairs of square Helmholtz coils cancel the Earth's magnetic field.
- Seven 37 mm 2nd-order gradiometers, 60 mm baseline.
- B_m ($\sim 200 \mu\text{T}$) continuously on; power supply with a large capacitor.
- Battery-powered current generators for gradients and spin-flip.
- B_p (45–65 mT) was battery-powered and ramped down through banks of solid-state switches.

Methods—Reference Channel




Ambient field deviation in the ULF MRI system location




- Manual ambient DC field compensation.
- A low-frequency dynamic cancellation system is being tested to enable automatic adjustments.
- Electronic compensation and software compensation have been tested and compared.
- In both cases it was possible to suppress noise lines from our NMR signals with a central Larmor frequency of 8.6 kHz.

Methods—Shielded System Hardware



- Wall-powered field generation
- 2nd-order gradiometers, 90 mm diameter & baseline
- Sensors co-exist with 100 mT B_p
 - Pb-Bi shields
 - 2nd-feedback
 - Compensation




Pb-Bi shield




Methods— B_p and Compensation




SQUIDs reach output limit w/o 2nd-feedback



With 2nd-feedback dynamic range maintained



Modeling—Parameters etc.




- Implemented in Matlab
- Bloch equations for NMR calculations
- Biot-Savart and reciprocity used for field calculations

	Tissue name	PD (%)	T_1 (ms)	T_2 (ms)
1	CSF	100	4360	329
2	GREY MATTER	86	635	83
3	WHITE MATTER	77	360	70
4	FAT	100	350	70
5	MUSCLE	100	120	47
6	MUSCLE/SKIN	100	120	47
7	SKULL	0	0	0
8	VESSELS	0	0	0
9	CONNECTIVE	77	500	61
10	DURA MATER	100	2569	329
11	BONE MARROW	77	500	70

Eleven high-resolution ($0.5 \times 0.5 \times 0.5 \text{ mm}^3$) volumes describing content of a voxel
 Montreal Neurological Institute, McGill University
 B Aubert-Broche, AC Evans, and DL Collins,
NeuroImage, **32**(1), 138–45, 2006.

Modeling—100 mT, 2 fT/ $\sqrt{\text{Hz}}$



Imaging parameters:

B_p : 100 mT

Polarization inversion time: 750 ms

Polarization time: 750 ms

Delay time: 10 ms

Encoding time: 35 ms

Acquisition time: 70 ms

N_y (phase): 103

N_z (phase): 41

Readout gradient, G_x : 7.0 Hz/mm


Phase gradient, G_y : 7.0 Hz/mm

Phase gradient, G_z : 3.0 Hz/mm

Voxel size: 2.0 × 2.0 × 4.8 mm³

Noise: 1.80 fT/ $\sqrt{\text{Hz}}$

Modeling—250 mT, 1 fT/ $\sqrt{\text{Hz}}$



Imaging parameters:

B_p : 250 mT

Polarization inversion time: 750 ms

Polarization time: 750 ms

Delay time: 10 ms

Encoding time: 35 ms

Acquisition time: 70 ms

N_y (phase): 103

N_z (phase): 41

Readout gradient, G_x : 7.0 Hz/mm

Phase gradient, G_y : 7.0 Hz/mm

Phase gradient, G_z : 3.0 Hz/mm



Voxel size: $2.0 \times 2.0 \times 4.8 \text{ mm}^3$

Noise: **0.90** fT/ $\sqrt{\text{Hz}}$

Modeling—Contrast (Inversion Recovery)


$t_p = 750 \text{ ms}$, $t_{inv} = 750 \text{ ms}$

$t_p = 2750 \text{ ms}$, $t_{inv} = 100 \text{ ms}$




Results—Unshielded

Unshielded MR-image
(4 averages; 65 mT)



Gelatin-agar phantom



Shielded MR-image
(1 average; 100 mT)



2D imaging sequence: $t_p = 1$ s, $t_{\text{enc}} = 100$ ms, $t_{\text{read}} = 200$ ms.

Phase encoding: 57 steps, $|G_z|_{\text{max}} = 1.62$ Hz/mm.

Frequency encoding (readout), G_x : 1.63 Hz/mm.

Resolution: $\sim 3 \times 3$ mm².

Spin-flip pulses: hard, 4 ms.

Gelatin-agar mixtures:

$T_2 \sim 120$ ms for the surrounding











$T_2 \sim 300$ ms for inclusions

Results—Unshielded cont.

Raw signal


PFA-compensation

Electronic compensation




Results—Inside MSR

4 s polarization (67 min.; 5 slices total)



Simulation





15 mm slices

Results—Inside MSR


3.5 s polarization (80 min.; 7 slices)

Simulation




Results—Inside MSR, contrast

4 sec polarization
and no inversion




Pulse-sequence





4 s polarization and
0.1 s inversion

0.5 s polarization and
0.5 s inversion

Outlook—Shorter TE (Imaging Sooner)




Outlook—Shorter TE (simulations)




- Signal from more tissues
- Signal in subsequent echoes

Outlook—New 250 mT B_p Coil







Field @ head center: ~250 mT
Field @ head edge: ~180 mT



Outlook—High B_p Issues



Data from the Nb-
gradiometers in the
unshielded system



Outlook—Heating of Ta-grads




Summary

- The shielded system works robustly with filtered amplifiers for all fields and gradients
- We need to increase B_p but not increase noise
- Improve duty-cycle to decrease imaging time
- Unshielded imaging with static Earth's field compensation and reference channel de-noising demonstrated
- Promising preliminary results from heating of Ta-gradiometers obtained

Acknowledgements

We greatly acknowledge the funding from Los Alamos National Laboratory's Laboratory Directed Research and Development (LDRD) office through grants 20130121DR and 20130624ER.

Thank you for your time and attention!



Robert H. Kraus Jr.
Michelle A. Espy
Per E. Magnelind
Petr L. Volegov

**ULTRA-LOW FIELD
NUCLEAR MAGNETIC
RESONANCE**

A NEW MRI REGIME

OXFORD